

PAW  WAVE™



Improve Mobility



Promote Bonding



Reduce Tension



Decrease Stress

You might have heard of our human products by PADO. We've taken our expertise in percussion & vibration massage and developed the first ever pet massager. Why? Because we love our Pets!


Just like you, we do everything we can to keep our pets healthy and happy. We've worked with everyone from engineers, to veterinarians, to the animals themselves all to bring the benefits of massage to our furry friends.

Want to meet our dream team of humans and happy pets?

Learn more about us at:
www.pawwave.pet

PAW WAVE BUZZ

VIBRATION THERAPY BRUSH MASSAGER



Using the BUZZ



For best results, charge the BUZZ before use.
(Approximately 180 mins for Full Charge)

Insert charging cable into the port on the side of your device.

Plug the charger into an outlet. The LED on the unit will turn on. The white LED lights will indicate charge level when unit is plugged into the charger. The LED indicator will turn off when the battery is fully charged.



Press and hold the power button to turn on the BUZZ.

Press the power button to adjust vibration mode speed. White LED lights on the massager will now indicate speed level.



Press and hold the power button to turn off the BUZZ.

Safety Feature: There is a 15 minute automatic shut-off for safety. If you need to use it longer, simply power it back on and massage away.



Calibrated for Pets

BUZZ was specifically designed for pets to provide an effective, yet gentle, therapeutic massage using vibration therapy. The custom silicone pad glides over fur and removes for easy cleaning. Vibration massage therapy can help relieve muscle tension around delicate tissue or sore spots. And it's a great way to help heal or relax your pet, as well as increase your bond with your best friend.

How to Attach

Line up Removable Therapy Brush Pad and cover the Brush of the massager.



Brush and Therapy Brush Pad

BUZZ was made for introducing massage to pets. The removable therapy brush pad is soft and easy to clean, offering a gentle yet soothing massage. When the pad is removed the hard plastic bristles offer a more penetrating massage, mimicking a deep and comfortable scratch.

WHAT'S IN THE BOX?



BUZZ Massanger



Removable Therapy Brush Pad



Charger & Cable



Pouch

For FAQs and other tips for using your
PAW WAVE PERK visit:

www.pawwave.pet

PADO INC.

28340 Ave Crocker, Suite #100
Santa Clarita, CA 91355